



ken-a-vision®
KNOWLEDGE THROUGH VISION



The Professor Stereo

ESH200 & ESH210

The worlds first inclined binocular stereo microscope designed specifically for students. This stereoscope is made of high impact resistant plastic and all glass optics for clean, crisp viewing. The Professor Stereo is available for purchase as a single unit or in packs of ten.

Additional Features and Benefits

- High Impact Plastic Construction
- Battery Operated - 300 hours of use
- Large Knobs for Small Hands
- Ken-A-Vision Cool Lighting
- Slip Clutch Indestructible Focus System
- Portable
- Complete with Dust Cover
- Student Proof with Safety Screw Construction

PROFESSOR STEREO SPECIFICATIONS

Weight	1.25 lbs / 0.57 kgs	Stage Plate	Reversible black / white
Size	5.5" X 4.5" X 9"	Upper Illuminator	Ken-A-Vision Cool Lighting
	14 cm X 12 cm X 23 cm	Base	5.5" x 4.5" / 14cm x 12cm
Eyepiece	10x	Power Source	2 AA Batteries (included)
Optics	Glass	ESH200	Single Unit
Head	45° angle fixed	ESH210	10 Stereo Microscopes
	Adjustable Interpupillary	Warranty	1 Year
Objectives	Binocular		
	2x		

Specifications subject to change without notice



ken-a-vision®
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com

kav.fly.ESH200.706v4