



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

Cordless CoreScope **T-1701C Series**

The Ken-A-Vision Cordless CoreScope has all the features and benefits of the Powered by Ken-A-Vision Cordless Microscopes. Plus, it has a light rheostat, a multi-directional round floating stage, and separate coarse and fine focusing knobs which make it easier for the student to operate.



Features and Benefits

- Light Rheostat
- Prism Head
- Monocular Head
- Rechargeable NiMH Batteries
- **Multi-Directional Floating Stage**
- Separate Coarse and Fine Focusing Knobs



T-1701C SPECIFICATIONS

Weight	4lbs / 2kgs	Stage	Multi directional floating stage with spring loaded stage clips
Size	12"Hx 6"L x 4"W	Light Source	Ken-A-Vision Cool Lighting
Prism Head	45° angle, rotates 360°	Light Board	Light Rheostat
Eyepiece	10x Widefield	Power Source	100,000 hours
Objective Lens	4x, 10x, 40xR DIN	Power Source Life	Rechargeable NiMH batteries
Stage Size	11.78" (30cm) Round	Warranty	500 cycles
Head	Monocular		(approximately 3 years)
Coaxial Focusing	Separate coarse & fine focusing knobs		10 Years

Specifications subject to change without notice



ken-a-vision[®]
KNOWLEDGE THROUGH VISION

(816) 353-4787 Fax (816) 358-5072

email: info@ken-a-vision.com

kav.fly.t1701C.0809v1