




**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

# Video Flex<sup>®</sup> 7100

Make learning fun with the Video Flex 7100! With a magnification power of up to 50x, the 7100 can even replace a microscope for many applications. The Ball and Socket head provides maximum flexibility. Theft protection is provided through a Kensington<sup>®</sup> security lock.



## Additional Features and Benefits

- ⚡ **Triangular Base** Slip and Tip proof
- ⚡⚡ **500 Lines of Resolution** For crisp clear images
- ⚡⚡ **Poly-Carbonate Build** Allows for rugged conditions and environment
- ⚡⚡ **C-mount Lens** Interchangeable lens system accommodates all c-mount lenses. 1/4 turn moves from 6mm to infinity
- ⚡⚡ **Quick Focus™** For recording or amplifying
- ⚡⚡ **Microphone** Adaptable to a wide range of lighting conditions
- ⚡⚡ **1.5 lux** Safety and Security
- ⚡⚡ **Kensington Slot**

### 7100 SPECIFICATIONS

Weight	5.95 lbs / 2.70 Kg	Resolution	500 Lines
Base Material	Polycarbonate	Magnification	50:1 and Greater
Base Dimension	7"x7"x7" / 18x18x18 cm	Exposure	Fixed Contrast
Height	30" / 76.2 cm	Output Signal	Composite Video
Head Type	Ball and Socket	Eyepiece Adapters	28mm & Built-in 34.5 mm
Lens	8 mm	Cable	12' A/V
Focus	¼ turn Quick Focus™	Warranty	5 Years
Focal distance	¼" / 6mm to Infinity		
Light Sensitivity	1.5 Lux		



**ken-a-vision**<sup>®</sup>  
KNOWLEDGE THROUGH VISION

(816) 353-4787 ⚡ Fax (816) 358-5072

email: [info@ken-a-vision.com](mailto:info@ken-a-vision.com)